


Tensoval®

comfort

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Preliminary remarks

 Please read these instructions carefully before first use as correct blood pressure measurement depends on the appropriate use of the device. These instructions for use are designed to instruct you, from the very start, in the individual steps of self-measurement of blood pressure using Tensoval comfort. You will thus receive important and helpful hints for producing reliable results for your personal blood pressure profile. Be sure to keep these instructions for use for future reference.

Table of Contents	Page
1. Introduction	132
2. General information on blood pressure	132
2.1 Significance of blood pressure values	132
2.2 Importance of self-measurement of blood pressure	133
2.3 Objectives of self-measurement of blood pressure	134
2.4 Regular blood pressure measurement	134
3. Getting ready for self-measurement	135
3.1 Inserting / changing the batteries	135
3.2 Setting date and time	135
3.3 10 golden rules for blood pressure measurement	135
3.4 Applying the cuff	137
4. Measuring blood pressure	138
5. Memory function	139
6. Explanation of error displays	140
7. Significance of symbols on the device and cuff	143
7.1 Control displays	143
7.2 Symbols	144
8. Important notes	144
8.1 Drugs	144
8.2 Pregnancy	144
8.3 Diabetes, history of other medical conditions	144
8.4 Arrhythmias, heart rhythm disorders, cardiac pacemakers	144
8.5 Important notes for self-measurement	145
9. Maintenance of the device	145
10. Accessories and spare parts	146
11. Warranty conditions	146
12. Technical data	147

	Page
13. Power supply, disposal notes and safety information	149
13.1 Batteries, mains adapters and disposal	149
13.2 Safety information	149
14. Legal requirements and guidelines	150
15. Instructions for the calibration check	150
16. Contact information for customer queries	151

1. Introduction

Congratulations on your purchase of this HARTMANN quality product. Tensoval comfort is a fully automatic blood pressure monitor for self-measurement on the upper arm which uses intelligent measuring technology (Fuzzy Logic) to take rapid and reliable measurement of systolic and diastolic blood pressure as well as pulse rate using an oscillometric measuring method.

2. General information on blood pressure

2.1 Significance of blood pressure values

To determine your blood pressure you need to measure two values:

- Systolic (upper) blood pressure: is produced when the heart contracts and pumps blood into the blood vessels.
- Diastolic (lower) blood pressure: this is the value measured when the heart muscle is dilated and again fills with blood.
- Blood pressure readings are expressed in mmHg.

The World Health Organisation (WHO) and the International Society of Hypertension (ISH) have developed the following classification for blood pressure values:

Assessment	Systolic pressure	Diastolic pressure
Optimal	up to 120 mmHg	up to 80 mmHg
Normal	up to 130 mmHg	up to 85 mmHg
Normal limit values	130 – 139 mmHg	85 – 89 mmHg
Grade 1 hypertension	140 – 159 mmHg	90 – 99 mmHg
Grade 2 hypertension	160 – 179 mmHg	100 – 109 mmHg
Grade 3 hypertension	over 180 mmHg	over 110 mmHg


Please note that this classification of blood pressure values is independent of age.

Established hypertension (high blood pressure) is defined as measurement of a systolic value greater than 140 mmHg and/or a diastolic value greater than 90 mmHg.

In general, blood pressure is considered to be too low in women (hypotension) with values of less than 100 mmHg systolic and less than 60 mmHg diastolic, and with values of less than 110 mmHg systolic and less than 70 mmHg diastolic in men. Please note that, unlike too-high blood pressure values, too-low blood pressure values are not usually expected to be associated with health risks.

2.2 Importance of self-measurement of blood pressure

Constantly elevated blood pressure multiplies the risk for other health problems. The most common causes of death worldwide are physical consequences such as heart attack, stroke and organic damages. Daily blood pressure monitoring is thus an important measure which will help to protect you from these risks.


2.3 Objectives of self-measurement of blood pressure

Your personal blood pressure profile is the regular documentation of the measured values over a certain period of time and therefore important information. In case of drug treatment in high blood pressure, your doctor can use your blood pressure profile to tailor your treatment to your particular needs. The better your drug treatment is tailored to your needs, the better you will feel. Regular, accurate blood pressure monitoring with Tensoval comfort will help you achieve this goal.

i Many people manage to lower their blood pressure through life-style changes (such as losing weight, dietary modification and getting more exercise) to levels that do not require drug treatment.

2.4 Regular blood pressure measurement

Numerous factors including physical exertion, taking drugs or the time of day may have an impact on blood pressure. Blood pressure should therefore always be measured at the same time of day under similar conditions.


i Our heart may beat up to 100,000 times a day, producing 100,000 different blood pressure values.

3. Getting ready for self-measurement

3.1 Inserting / changing the batteries

Open the battery cover on the underside of the device (see fig.) Insert batteries (see chap. 12, "Technical data"), ensuring correct polarity ("+" and "-"). Close the battery lid.

3.2 Setting date and time

- After inserting the batteries the time function automatically appears. The flashing year is then displayed. The year can be changed using the M1 (+) and M2 (-) buttons. Store the year by pressing the START/STOP button.
- Next the month is stored. The number on the right will flash. Use the buttons as above for setting the year. Repeat procedure for saving the day, hour and minute. Each time you change the batteries the date and time must be reset.

3.3 10 golden rules for blood pressure measurement

Many factors are involved when measuring blood pressure. These ten general rules will help you to take the readings correctly.

1. Rest for approx. 5 minutes before measurement. Even deskwork increases blood pressure by an average of approx. 6 mmHg systolic and 5 mmHg diastolic.

2. Do not consume any coffee or nicotine up to one hour before measurement.

3. Do not measure when you have a strong urge to urinate. A full bladder can lead to an increase in blood pressure of approx. 10 mmHg.

4. Take measurements from the naked upper arm and while sitting upright.

5. In the case of using a wrist monitor, hold the cuff at the level of the heart


during the measuring procedure. The cuff of an upper arm monitor automatically finds the correct level at the arm.

6. Do not talk or move during the measuring procedure. Talking increases the values by approx. 6 – 7 mmHg.


7. Wait at least one minute between two measurements, so that the vessels are relieved from pressure in preparation for a new measurement.

8. Enter values in the blood pressure diary: note down the measured values, together with any drugs taken, the date and time in your blood pressure diary.

9. Take measurements regularly. Even if your values have improved, you should continue to check them for monitoring purposes.


10. Always take measurements at the same time of day. Because a person has approx. 100,000 different blood pressure values every day, individual measurements have no significance. Only regular measurements at the same time each day over a long period of time allow a meaningful evaluation of blood pressure values.

**Further notes:**

- You should take your blood pressure in a quiet place, in a relaxed seated position. Measurement can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
- Do not take your blood pressure after taking a bath or exercising.

the lower edge of the cuff is approx. 2.5 cm from the bend of the elbow (see fig.).

**3.4 Applying the cuff**


- Measurement should be taken on the naked arm which usually has the higher blood pressure value.
- When using the draw-clamp cuff (standard cuff), included with the device, thread the end of the cuff through the metal hoop, so as to form a loop. In this instance, the Velcro fastener must be on the outside. Wrap the cuff round the upper arm, the air tube lies in the centre of the elbow angle, running down on the inside of the forearm on a level with the middle finger, pointing towards the hand. Take the free end of the cuff, wrap it firmly round the arm and close the Velcro fastener.
- Check that the white strip is placed in the middle of the elbow, on the artery and that
- The cuff should be firm but not too tight. You should be able to push two fingers between the arm and the cuff. Please note that uneven wrapping of the cuff can lead to inaccurate readings.
- Use the markings on the edge of the cuff to check whether cuff size is correct. The white arrow should lie within the red marking strip.
- Insert the cuff connector in the cuff socket on the left side of the device. Attention: do not use the mains adapter socket on the back of the device!
- Make sure that the tube is not kinked or damaged, as the measurement could otherwise be impaired.

4. Measuring blood pressure

- We recommend that you measure your blood pressure while sitting with your back supported by the back of the chair. Place both feet flat on the floor with the legs next to each other. Place your forearm with the palm relaxed upwards on a support and make sure the cuff is at the level of the heart.
- Do not turn the device on until the cuff has been applied, otherwise the cuff can become damaged through the resulting excess pressure.
- Press the START/STOP button. The appearance of all display segments followed by a flashing arrow pointing downwards, shows that the device is checking itself automatically and is ready for use. Subsequently, automatic inflation begins at approx. 180 mmHg. If the inflation pressure is insufficient or if the measurement is interrupted, the device continues to inflate at a rate of 30 mmHg until a high enough pressure is reached.
- If you require a higher inflation pressure, you can avoid having to repeat inflating by pressing the START/STOP button again shortly after inflation starts for a few seconds and keep it pressed until the desired cuff pressure is reached. This should be approx.

30 mmHg over the systolic value.

- Important: you should not move or talk throughout the entire measuring procedure!
- As pressure in the cuff decreases, the heart symbol and the falling cuff pressure are displayed.
- A beep indicates the end of measurement. Then the systolic and diastolic blood pressure values appear simultaneously on the display, with the pulse rate beneath them (see fig.).

The time appears above the reading and M1 or M2 is displayed on the left. M1 represents the reading for the first person and M2 for the second. As long as the reading is displayed, you can apply the values to the respective person by pressing the M1 or M2 button. Otherwise, if you do not apply them the reading is automatically stored for the memory button appearing on the display.

- In order to switch off the device, press the START/STOP button. Otherwise the device will switch

itself off automatically after 3 minutes.

- If you wish to stop measurement for any reason, simply press the START/STOP button. The inflation or measuring procedure is interrupted and an automatic fall in pressure occurs.
- If this symbol  can be seen to the left of the display next to the diastolic reading (DIA), the device has detected an irregular heart-beat during measurement. However, the measurement may also have been disrupted by body movement or speaking. It is best to repeat the measurement. If you regularly see this symbol during your blood pressure measurements, we recommend having your heart rhythm checked by your doctor.

5. Memory function

- The memory recall is activated by pressing the memory button, when the device is switched off. Press M1 for the first person's stored values and M2 for the second person. The corresponding symbol M1 or M2 will be displayed. First, the average value of all stored data will be shown for the corresponding person. An "A" appears on the display and the number of readings, from which the average value was calculated, is displayed in the top right-hand corner (see fig.).

Press the memory button to access memory position number 1.

- Tensoval comfort can store up to 60 measurements per memory button. The most recent measured value is always at memory position number 1 and all old measured values shift down one memory position. When all memory positions are occupied the oldest value will be deleted each time.
- By repeatedly pressing the memory button all stored values can be recalled one by one.
- The measured value and the relevant number of the memory position appear when you access stored values. The memory position, date and time appear in 2 – 3 second intervals.
- You can cancel the memory function at any time by pressing the START/STOP button. Otherwise the device will switch itself off automatically after a few seconds.
- Even if the power supply fails, e.g. when changing the batteries, the stored values will still be available.

i If an irregular heartbeat was detected during measurement, this information  is also stored and displayed when recalling the measured values from the device memory, together with the systolic and diastolic blood pressure value, the pulse rate, the time, date and year.

Deleting the stored values

You can delete all stored data in M1 or in M2 for the respective

person separately. To do this, press the memory button of the corresponding person. The average value will then appear on the display. Press down on the memory button again and hold it down for a longer time. After four seconds the display will flash and after eight seconds all data for the corresponding person will be deleted and only M1 or M2 will be displayed. If you release the memory button ahead of time, no data will be deleted.

6. Explanation of error displays

Error which has occurred	Possible causes	Remedy
Device will not turn on.	No batteries, they have been inserted incorrectly or are dead.	Check batteries and insert four identical, new batteries, if and when necessary.
	Mains adapter not correctly connected or defective.	Ensure the mains adapter is plugged in the connecting socket on the rear side of the device.
Cuff will not inflate.	Cuff connector is incorrectly positioned in the connecting socket of the device.	Check connection between the cuff connector and the connecting socket.
	Wrong cuff type connected.	Check to make sure that only the approved Tensoval comfort cuffs and the corresponding connectors were used.

Error which has occurred	Possible causes	Remedy
	<p>The measuring signals could not, or not correctly, be read. This can be caused by incorrect application of the cuff, moving, talking or by a very weak pulse.</p>	<p>Check the correct position of the cuff. Do not talk or move during the measuring procedure. Also observe the 10 golden rules as well as the notes in the box in chapter 3.3.</p>
	<p>Cuff does not inflate or does not inflate fast enough. This can, among other things, be due to a too loosely applied cuff or to movement.</p>	<p>Apply the cuff so that you are able to push about two fingers between the cuff and the upper arm. Air tube is not correctly inserted into the device.</p> <p>Check the correct position of the connector. If this error occurs often, use a new cuff.</p>
	<p>Air release during the measuring procedure is too fast or too slow. The cuff could have become undone or loosened itself. A movement during the measuring procedure is also a possibility.</p>	<p>Check the correct position of the cuff. Do not move during the measuring procedure.</p>
	<p>The pressure in cuff exceeds 300 mmHg. Therefore, an automatic fall in pressure occurs.</p>	<p>Please rest for at least a minute and take the measurement again.</p>

Error which has occurred	Possible causes	Remedy
	If the battery symbol is flashing the batteries are almost dead. Only a few more measurements are possible.	Keep new batteries of the same type handy (type AA LR06).
	If the battery symbol is permanently illuminated the batteries are dead and have to be replaced.	Insert new batteries of the same type (type AA LR06). However, please note that the battery sign which appears when the device is turned on together with all display functions, and then goes out again, does not provide an indication of the battery charging level.
Implausible measured values.	Incorrect cuff size.	Use the cuff corresponding to your upper arm size.
	Cuff placed on top of clothing.	Apply cuff on the naked skin.
	Rolled-up clothing impedes blood circulation.	Wear loose clothing. Make sure that rolled-up sleeves do not impair circulation in the upper arm.
	Cuff wrongly applied.	Take note of the instructions and images showing how to apply the cuff correctly to the upper arm.
	Cuff tube folded or squashed.	Ensure that the cuff tube lies straight and loose.
	Cuff was not correctly inflated.	Check the correct position of the upper arm cuff.

Error which has occurred	Possible causes	Remedy
	Moving, talking or excitement during the measuring procedure.	Please take measurements in a relaxed position whilst seated. Do not talk or move during the measuring procedure.
	Lack of relaxation before taking a measurement.	Relax for 5 minutes before taking a measurement.
	Stimulants taken before measurement.	Please avoid alcohol / nicotine and caffeine for one hour before taking a measurement.

Switch the device off if an error symbol appears. Check the possible causes and note the 10 golden rules in chapter 3 and the instructions in chapter 8. Relax for a minute and take the measurement again. Do not move or talk during the measurement.

7. Significance of symbols on the device and cuff

Appears during inflation.

Appears during automatic checking.

7.1 Control displays

Flashes when the device is measuring and the pulse is being taken.

Displays the stored measured values for person 1.

Replace batteries if the battery symbol is permanently illuminated.

Displays the stored measured values for person 2.

Measuring error, cf. chap. 6.

7.2 Symbols

Pay attention to the operating instructions.

Please note.

Protection from electric shock (type BF).

Disposal note.

8. Important notes

8.1 Drugs

Self-measurement of blood pressure does not replace treatment! So do not interpret your measured values on your own and do not use them for self-prescribed treatment. Take measurements as instructed by your doctor and have confidence in his diagnosis. Take drugs as prescribed by your doctor and never alter the dose on your own. Discuss the appropriate time for self-measurement of blood pressure with your doctor.

8.2 Pregnancy

Blood pressure may change during pregnancy. Regular blood pressure monitoring is particularly important if you have high blood pressure because the elevated blood pressure values may affect the development of the foetus. Check with your doctor whether and, if so, when you should carry out self-measurement of blood pressure.

8.3 Diabetes, history of other medical conditions

If you have diabetes, hepatic disorders or narrowed blood vessels (e.g. arteriosclerosis, peripheral arterial occlusive diseases (PAOD), you should consult your doctor before carrying out self-measurement because altered measured values may occur in such cases. If you suffer from certain blood diseases (e.g. haemophilia), severely impaired blood flow, or if you take blood-thinning drugs, you should also ask your doctor before carrying out self-measurement.

8.4 Arrhythmias, heart rhythm disorders, cardiac pacemakers

■ In the case of severe heart rhythm disorders (arrhythmias) measurements should only be taken in consultation with the doctor. Due to the oscillometric measuring method, in some cases incorrect readings may be

determined or no measurement results are obtained (Err).

- If this symbol  appears often, it can be a sign of heart rhythm disorders. In this case, consult your doctor. Severe heart rhythm disorders may produce false measurements or impair the measuring accuracy. Please discuss with your doctor whether self-measurement of blood pressure is suitable for you. Altered measured values may occur in the case of self-measurement performed by cardiac pacemaker wearers. The blood pressure monitor itself has no impact on the cardiac pacemaker. Please note that the displayed pulse rate is not suitable for checking the rate of cardiac pacemakers. Please check with your doctor whether self-measurement of blood pressure is advisable if you are wearing a cardiac pacemaker.

8.5 Important notes for self-measurement

- Individual readings are situation-related and thus are not useful.
- Even slight changes in internal and external factors (e.g. deep breathing, stimulants, talking, excitement, climatic factors) lead to fluctuations in blood pressure. This is why your doctor and pharmacist often obtain different readings.

- Always measure blood pressure on the same arm and rest the forearm relaxed on a support.
- Measurement can be taken on the right or left arm. The arm giving higher readings should be used for long-term blood pressure monitoring.
- Patients who tend to bruise and/or are sensitive to pain on pressure, should only take blood pressure measurements after consulting a doctor.
- The medium cuff is suitable for an upper arm circumference between 22 and 32 cm, the large cuff is suitable for an upper arm circumference between 32 and 42 cm. Outside these limits correct readings may no longer be guaranteed.

9. Maintenance of the device

- Do not expose the device neither to extreme temperatures nor to humidity, dust, or direct sunlight because this may lead to malfunction.
- This device consists of high-quality electronic precision components. Protect the device from knocks and do not immerse in water.
- The moulded cuff (accessory) should not be folded or over-stretched.

- Never open the device. Repairs may only be carried out by authorized professionals.
- Only use a soft, moistened cloth to clean the device. Do not use detergents or solvents.
- The cuff can be cleaned carefully with a lightly moistened cloth and mild soap solution. Do not completely immerse the cuff in water.

10. Accessories and spare parts

To ensure measurement accuracy, only use original HARTMANN accessories which may be obtained from your pharmacist or specialist medical supplier.

Upper arm circumference	Required cuff
22 – 32 cm	medium
32 – 42 cm	large

Moulded cuff, medium for upper arm circumference of
22 – 32 cm
Code no. 900 166

Draw-clamp cuff (standard cuff), medium for upper arm circumference of
22 – 32 cm
Code no. 900 154

Draw-clamp cuff (standard cuff), large for upper arm circumference of
32 – 42 cm
Code no. 900 155

Tensoval mains adapter
Code no. 900 152

11. Warranty conditions

We give a 3-year warranty on this high-quality device for measuring blood pressure from the day of purchase and in accordance with the following conditions.

Claims must be made during the warranty period. The date of purchase may be documented by the appropriately completed and stamped warranty document or proof of purchase.

Within the warranty period, HARTMANN shall replace or repair any faulty device components free of charge which were caused by material or manufacturing errors. This does not extend the warranty period.

This warranty is not applicable to damage caused by improper use or unauthorized interference. Parts that are subject to wear and tear (batteries, cuffs, mains adapters etc.) are excluded from the warranty. Claims for compensation

are limited to the value of the goods; compensation for subsequent damages is expressly excluded.

In warranty cases please send the device with cuff and, if applicable, the mains adapter together with the fully completed and stamped warranty certificate direct, or via your dealer to the Customer Services department for your country.

AE – PAUL HARTMANN
Middle East FZE
Dubai

AU – PAUL HARTMANN Pty. Ltd.
Level 6, 5 Ryder Boulevard
Rhodes, NSW 2138 Australia

HK – PAUL HARTMANN
Asia-Pacific Ltd.
Hong Kong

ZA – HARTMANN South Africa
2194 Johannesburg

12. Technical data

Measuring method: oscillometric

Display range: 0 – 300 mmHg

Measuring range: Systole (SYS): 50 – 250 mmHg
Diastole (DIA): 40 – 180 mmHg
Pulse: 40 – 160 beats / minute

The displaying of values outside the measuring range cannot be guaranteed.

Technical measuring accuracy:

Cuff pressure: +/- 3 mmHg
Pulse: +/- 5% of displayed pulse rate

Clinical measuring accuracy:

complies with the European Standard EN 1060, Part 3

Power supply:

4 x 1.5 V Mignon alkaline-manganese (AA/LR06) batteries or optional HARTMANN mains adapter

Battery capacity:

Tensoval comfort: > 1,500 measurements
Tensoval comfort large: > 1,000 measurements

Protection from electric shock:	Protection class II (when using the Tensoval mains adapter) Medical electrical equipment with an internal power supply (when using batteries) Applied part: type BF
Protection against harmful penetration of water or solid materials:	IP20
Operating mode:	continuous operation
Inflation pressure:	approx. 180 mmHg
Automatic switch-off function:	3 minutes after end of measurement
Cuff:	Standard cuff, medium 22 – 32 cm Standard cuff, large 32 – 42 cm Moulded cuff, medium (optionally) 22 – 32 cm
Pressure release valve:	electronically controlled linear valve
Memory capacity:	2 x 60 measurements and mean value
Operating conditions:	Ambient temperature: + 10 °C to + 40 °C (+ 50 °F to + 104 °F) Relative humidity: 15 – 85 %
Storage / transport conditions:	Ambient temperature: - 20 °C to + 50 °C (- 4 °F to + 122 °F) Relative humidity: 15 – 85 %
Serial number:	In battery compartment
	PAUL HARTMANN AG Paul-Hartmann-Straße 89522 Heidenheim/Germany
	 0 1 2 3

13. Power supply, disposal notes and safety information

13.1 Batteries, mains adapters and disposal

- We recommend the use of high-quality batteries as other batteries or accumulators may result in a reduction in the measuring performance. Never mix old and new batteries or batteries made by different manufacturers.
- If the battery symbol is permanently displayed you should change the batteries as soon as possible. Please note that the symbol always looks “empty”.
- Remove the batteries from the device if it is not being used for a longer period.
- In the interests of environmental protection exhausted batteries may not be disposed of in household waste. Please observe the applicable waste disposal regulations or use public collecting bins. 
- Operation with mains adapter: on the rear of the device there is a connection socket for the mains adapter (output voltage 6V DC/600mA). Please use only a HARTMANN mains adapter, which may be obtained from your pharmacist or specialist medical supplier. Otherwise the

measuring accuracy of the device cannot be guaranteed.

13.2 Safety information

- The device is not waterproof!
- Do not leave the device unattended near toddlers or persons who cannot operate it themselves.
- Use the device for taking blood pressure measurements on the upper arm only.
- Do not under any circumstances carry out blood pressure measurements on babies or toddlers.
- Do not expose the device to hard knocks or vibrations.
- Do not drop the device to the floor.
- Do not excessively bend or fold the arm cuff and the air tube.
- The device must not be altered, dismantled, or repaired by the user.
- Use the device only with the approved upper arm cuff; otherwise the device can be damaged externally or internally.
- The cuff tube may only be removed from the device by pulling the black connector. Never pull on the tube itself.
- Never inflate the cuff when it is not properly applied to the upper arm
- Please do not apply the cuff over a wound, as this may result in further injuries.

- If you have had a mastectomy, do not carry out the measurement on the arm on the affected side of the body.
- Please note that the pressure built up by the cuff can lead to temporary disruption to medical devices being simultaneously used on the same arm.
- If an intravenous treatment is being carried out or a venous catheter is present on the arm, blood pressure measurements can lead to injury. Never use the cuff on the arm on which these conditions apply.
- Please wait for one minute between two measurements and ensure that the tube is not knotted, kinked, or damaged.
- If you are carrying out the measurement on another person, please ensure that the use of the device does not result in persistent impairment of the blood circulation.

14. Legal requirements and guidelines

Tensoval comfort complies with the requirements of the EC directive 93/42/EEC on medical devices (Medical Device Directive MDD) and bears the CE mark.

The device complies, for example, with the European Standard EN 1060:

Non-invasive blood pressure measuring devices, Part 1: General requirements and Part 3: Additional requirements for electro-mechanical blood pressure measuring systems.

Clinical testing of measuring accuracy was performed according to the European Standard EN 1060-4. The requirements according to the ANSI/AAMI test protocol SP10-1992 are also fulfilled.

Portable and mobile high-frequency and communication devices, such as telephone and mobile phone, can impair the functional capability of electronic medical devices. In compliance with European Standard EN 60601-1-2 further information can be supplied by HARTMANN on request.

Over and beyond the legal requirements, the device has been validated by the ESH (European Society of Hypertension) in accordance with the ESH-IP2 protocol.

15. Instructions for the calibration check

We recommend a calibration check at intervals of two years in the case of professionally used devices, for example in pharmacies, medical practices, or hospitals. You should also observe the national regulations determined by the legislator,

such as, in Germany, the “Medizinprodukte-Betreiberverordnung” (Medical Device Operating Regulation). Calibration checks can be carried out either by competent authorities or authorised maintenance providers against compensation.

Instructions for the calibration mode:

Remove the batteries in order to switch to calibration mode. Press down the START/STOP button and then insert the batteries again. Release the button after a few seconds and, after a few moments, two zeros will appear one above each other on the display. Instructions on the calibration check will be supplied on request to competent authorities or authorised maintenance providers by HARTMANN.

16. Contact information for customer queries

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
■ English

- 1 Time / date
- 2 Systolic value
- 3 Diastolic value
- 4 Pulse
- 5 Battery symbol
- 6 Appears during automatic checking
- 7 Appears during inflation
- 8 Flashes when the device is measuring and the pulse is being taken
- 9 Irregular heartbeat
- 10 Memory user 2
- 11 Memory user 1

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PAUL HARTMANN AG · 89522 Heidenheim, Germany