HARTMANN

Cosmos® hydro

Hydrocolloid plasters for blisters

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EN Instructions for use

Product description

The hydrocolloid plaster is a non-sterile, single-use, self-adhesive plaster with a 100 % hydrocolloid* formulation and occlusive properties (moist treatment). Plasters are cushioning, waterproof.

Composition

The plaster is composed of two layers of adhesive hydrocolloid material (containing carboxymethyl cellulose): one acting as functional pad and the other acting as the dressing backing material. *Functional pad and backing material are made of 100 % hydrocolloid. Silicone paper protects the side of the plaster with adhesive coating.

Properties

The main function of the product is the provision and maintenance of optimal hydration of the exposed tissues, by the moist wound healing properties of the hydrocolloid formulation, and the protection of blisters, by the adhesion properties of the backing material. The self-adhesive hydrocolloid dressing contains a layer which absorbs exudate and expands into a gel-like mass, thus providing cushioning for the treated damaged area and helping to relieve pressure and friction.

Intended use

Hydrocolloid plasters are non-active, single-use, self-adhesive, non-sterile devices consisting of a self-adhesive hydrocolloid wound contact layer and an adhesive coated backing material (carrier material), protected with a silicone paper, for the protection and treatment of blisters and absorption of exudates. The products are intended for humans, without restrictions (in age group, weight range, health status or condition), for use by lay users, in contact with intact skin (healthy skin) and breached skin. The products do not contain resorbable components, medicinal products/ drug substances (derived from human blood or plasma or not), tissues or cells of human or animal origin, radioactive substances, nanomaterials or other hazardous components.

Purpose

Cosmos® hydro hydrocolloid plasters for blisters are designed to effectively protect blisters from contamination, wetting, exposure to the external environment, bacteria, maintain a natural level of moisture, absorb exudate from calluses, reduce pain, create optimal conditions for healing, and also to prevent the formation of blisters.

Indication for use

Protection of blisters from contamination, wetting, exposure to the external environment, and bacteria.
Treatment of blisters.
Maintaining a natural level of moisture.
Absorption of exudate.

Contraindications for use

- untreated infected blisters,
- untreated infected minor (superficial) skin lesions

- serious injuries, severely bleeding or profusely exuding wounds, bedsores, or for any purpose other than those listed in the intended use

- individual intolerance to manufacturing materials.

Possible side effects

Allergic reaction to the materials used to make the product.

Product disposal

The disposable components and used plasters can be disposed of normally with household waste.

Mode of application

- Wash hands and ensure the affected area is clean, apply the plaster only to dry and oil-free skin.
- Open the package and pull out the plaster.
- Warm the plaster with your hands before use.
- Remove the first protective paper and place the exposed part of the plaster to the affected skin area.
- Then gently pull the second protective paper and apply the rest of the plaster. Gently mould the plaster edges to ensure optimum adhesion.
- Products containing hydrocolloid can be applied for a maximum of 4 days for a single product. To remove the plaster early or after the maximum application time has been reached, moisten it briefly with warm water.

Special precautions

- For single use only. Not reusable.
- Apply the plaster only if the hydrocolloid pad covers the entire surface of the blister.
- Consult a physician before applying the plaster on a potentially infected blister or because signs of infection (eg, redness, heat, swelling, or pain) may appear.
- In case of diabetes mellitus, always consult your doctor about the use of the plaster.
- Do not puncture an unruptured blister and put yourself at risk of infection.
- Do not use the product for more than 4 days.
- Store the plasters in a cool (not refrigerated) and dry place at a temperature of +5 °C to +30 °C, keep out of the reach of children.

Incident reporting

For a patient/user/third party in the European Union and in countries with identical regulatory regime (Regulation (EU) 2017/745 on Medical Devices); if, during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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